

Anti-equine IL-1 alpha Polyclonal Antibody

KPAB0131

Description

This Anti-equine IL-1 alpha Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0131
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Goat
Reactivity:	Equine
Immunogen:	Recombinant Horse IL-1 alpha
Molecular Weight:	18.1 kDa (calculated) for equine IL-1 α
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Bovine IL-1A: Moderate; Canine IL-1A: Strong; Feline IL-1A: Strong
Homology:	Equus caballus (horse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.