

Biotinylated Anti-chicken IL-21 Polyclonal Antibody

KPAB0141

Description

This Biotinylated Anti-chicken IL-21 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0141
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Goat
Reactivity:	Chicken
Immunogen:	Recombinant Chicken IL-21
Molecular Weight:	14.2 kDa (calculated) for chicken IL-21
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Homology:	Gallus gallus (chicken)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.