

Anti-dolphin IFN gamma Polyclonal Antibody

KPAB0145

Description

This Anti-dolphin IFN gamma Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0145
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Goat
Reactivity:	Dolphin
Immunogen:	Recombinant Dolphin IFN gamma
Molecular Weight:	17.2 kDa (calculated) for dolphin IFN- η
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Canine IFN-H: Weak; Equine IFN-H: Weak; Human IFN-H: Moderate; Mouse IFN-H: Weak
Homology:	Orcinus orca (killer whale), Tursiops truncatus (bottlenose dolphin)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.