

Anti-bovine IL-6 Polyclonal Antibody

KPAB0153

Description

This Anti-bovine IL-6 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0153
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Chicken
Reactivity:	Bovine
Immunogen:	Recombinant Bovine IL-6
Molecular Weight:	20.7 kDa (calculated) for bovine IL-6
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN3
Cross-Reactivity:	Caprine IL-6: Strong; Ovine IL-6: Strong; Rat IL-6: Weak
Homology:	Bos indicus (zebu), Bos javanicus (banteng), Bos taurus (cattle)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.