

## Biotinylated Anti-bovine IL-6 Polyclonal Antibody

### KPAB0154

#### Description

---

This Biotinylated Anti-bovine IL-6 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

#### Product Information

---

<b>SKU:</b>	KPAB0154
<b>Concentration:</b>	1.0 mg/mL
<b>Conjugated:</b>	Biotin
<b>Applications:</b>	ELISA, ELISpot

#### Antibody Data

---

<b>Host:</b>	Chicken
<b>Reactivity:</b>	Bovine
<b>Immunogen:</b>	Recombinant Bovine IL-6
<b>Molecular Weight:</b>	20.7 kDa (calculated) for bovine IL-6
<b>Purification:</b>	Antigen-affinity Purification
<b>Formulation:</b>	1X PBS buffer containing 0.09% NaN3
<b>Cross-Reactivity:</b>	Caprine IL-6: Strong; Ovine IL-6: Strong; Rat IL-6: Weak
<b>Homology:</b>	Bos indicus (zebu), Bos javanicus (banteng), Bos taurus (cattle)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.