

Biotinylated Anti-canine IFN alpha Polyclonal Antibody

KPAB0166

Description

This Biotinylated Anti-canine IFN alpha Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0166
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Rabbit
Reactivity:	Canine
Immunogen:	Recombinant Canine IFN alpha
Purification:	Antigen-affinity Purification
Homology:	Canis lupus dingo (dingo), Canis lupus familiaris (dog)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.