

Biotinylated Anti-equine IFN alpha 1 Polyclonal Antibody

KPAB0168

Description

This Biotinylated Anti-equine IFN alpha 1 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0168
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Rabbit
Reactivity:	Canine
Immunogen:	Recombinant Canine IFN alpha
Molecular Weight:	18.5 kDa (calculated) for equine IFN α 1
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Bovine IFN α a: Weak; Canine IFN-A: Weak; Swine IFN α 1: Moderate
Homology:	Equus caballus (horse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.