

## Biotinylated Anti-chicken IFN alpha Polyclonal Antibody

### KPAB0170

#### Description

---

This Biotinylated Anti-chicken IFN alpha Polyclonal Antibody has been validated for ELISA, ELISpot applications.

#### Product Information

---

<b>SKU:</b>	KPAB0170
<b>Concentration:</b>	1.0 mg/mL
<b>Conjugated:</b>	Biotin
<b>Applications:</b>	ELISA, ELISpot

#### Antibody Data

---

<b>Host:</b>	Goat
<b>Reactivity:</b>	Chicken
<b>Immunogen:</b>	Recombinant Chicken IFN alpha
<b>Molecular Weight:</b>	19.0 kDa (calculated) for chicken IFN- $\alpha$
<b>Purification:</b>	Antigen-affinity Purification
<b>Formulation:</b>	1X PBS buffer containing 0.09% NaN <sub>3</sub>
<b>Homology:</b>	Gallus gallus (chicken)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.