

Biotinylated Anti-equine IL-1 beta Polyclonal Antibody

KPAB0196

Description

This Biotinylated Anti-equine IL-1 beta Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0196
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Rabbit
Reactivity:	Equine
Immunogen:	Recombinant Horse IL-1 beta
Molecular Weight:	17.3 kDa (calculated) for equine IL-1 β
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Homology:	Equus asinus (ass), Equus caballus (horse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.