

Anti-bovine LIF Polyclonal Antibody

KPAB0199

Description

This Anti-bovine LIF Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0199
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Chicken
Reactivity:	Bovine
Immunogen:	Recombinant Bovine LIF
Molecular Weight:	19.8 kDa (calculated) for bovine LIF
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.02% NaN ₃
Cross-Reactivity:	Canine LIF: Weak; Rabbit LIF: Moderate; Swine LIF: Weak
Homology:	Bos taurus (cattle)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.