

Anti-equine CCL2 Polyclonal Antibody

KPAB0203

Description

This Anti-equine CCL2 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0203
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Goat
Reactivity:	Equine
Immunogen:	Recombinant Horse CCL2
Molecular Weight:	8.5 kDa (calculated) for equine CCL2
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Bovine CCL2: Weak; Canine CCL2: Moderate; Swine CCL2: Weak
Homology:	Equus caballus (horse), Equus przewalskii (Przewalski's horse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.