

Anti-rabbit SCF Polyclonal Antibody

KPAB0296

Description

This Anti-rabbit SCF Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0296
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Chicken
Reactivity:	Rabbit
Immunogen:	Recombinant Rabbit SCF
Molecular Weight:	20.7 kDa (calculated) for rabbit SCF
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Homology:	Oryctolagus cuniculus (rabbit)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.