

Anti-canine PD-L1 Polyclonal Antibody

KPAB0319

Description

This Anti-canine PD-L1 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0319
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Canine
Immunogen:	Recombinant Canine PD-L1
Molecular Weight:	24.7 kDa (calculated) for canine PD-L1
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Homology:	Canis lupus dingo (dingo), Canis lupus familiaris (dog)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.