

Biotinylated Anti-mouse CCL2 Polyclonal Antibody

KPAB0360

Description

This Biotinylated Anti-mouse CCL2 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0360
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit
Reactivity: Mouse
Immunogen: Recombinant Mouse CCL2
Purification: Antigen-affinity Purification
Homology: Mus musculus (house mouse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.