

Anti-swine SCF Polyclonal Antibody

KPAB0383

Description

This Anti-swine SCF Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0383
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Swine
Immunogen:	Recombinant Porcine SCF
Molecular Weight:	18.7 kDa (calculated) for swine SCF
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Bovine SCF: Weak
Homology:	Sus scrofa (pig)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.