

Anti-cynomolgus monkey CCL4 Polyclonal Antibody

KPAB0409

Description

This Anti-cynomolgus monkey CCL4 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0409
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Cynomolgus Monkey
Immunogen:	Recombinant Cynomolgus Monkey CCL4
Molecular Weight:	7.8 kDa (calculated) for cynomolgus monkey CCL4
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Bovine CCL4: Moderate; Canine CCL4: Moderate; Equine CCL4: Weak; Human CCL4: Strong; Swine CCL4: Weak
Homology:	Macaca fascicularis (crab-eating macaque), Macaca mulatta (Rhesus monkey), Rhinopithecus roxellana (golden snub-nosed monkey)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.