

Biotinylated Anti-ovine IL-4 Polyclonal Antibody

KPAB0422

Description

This Biotinylated Anti-ovine IL-4 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0422
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit
Reactivity: Ovine
Immunogen: Recombinant Sheep IL-4
Purification: Antigen-affinity Purification
Homology: Budorcas taxicolor (taken), Ovis aries (sheep)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.