

Biotinylated Anti-ovine IFN gamma Polyclonal Antibody

KPAB0430

Description

This Biotinylated Anti-ovine IFN gamma Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0430
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Rabbit
Reactivity:	Ovine
Immunogen:	Recombinant Sheep IFN gamma
Molecular Weight:	16.9 kDa (calculated) for ovine IFN- η
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Bovine IFN-H: Moderate; Caprine IFN-H: Moderate; Dolphin IFN-H: Weak
Homology:	Ovis aries (sheep)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.