

Glu-Glu Tag Monoclonal Antibody

MACO0426

Description

This Glu-Glu Tag Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	MACO0426
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Synonyms:	-
Clone:	-
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	All

Antibody Data

Isotype:	-
Uniprot:	N A
Host Species:	Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Glu-Glu Tag
Immunogen Species:	-
Buffer:	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:500-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-