

ZBTB45 Antibody

PACO06996

Description

This ZBTB45 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO06996
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	ZBTB45 antibody, ZNF499 antibody, Zinc finger and BTB domain-containing protein 45 antibody, Zinc finger protein 499 antibody
Clone:	-
Applications:	WB IHC
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse, Rat

Antibody Data

Isotype:	IgG
Uniprot:	Q96K62
Host Species:	Rabbit
Purification:	-
Immunogen:	Recombinant Protein
Immunogen Species:	Homo sapiens (Human)
Buffer:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:500-1000
IHC	1:200-500	

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol


Validation Data

Image



Description

Immunohistochemical analysis of paraffin-embedded Mouse Kidney Tissue using ZBTB45 Polyclonal Antibody.



Western blot analysis of 1) HeLa, 2) Jurkat, 3) Mouse Brain, 4) Mouse Kidney, 5) Rat Brain using ZBTB45 Polyclonal Antibody.