

## BEST3 Antibody

PACO08005

### Description

---

This BEST3 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	PACO08005
<b>Contents:</b>	100µl Bradford Reagent: 1 vial (2ml)
<b>Category:</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Synonyms:</b>	BEST3, VMD2L3, Bestrophin-3, Vitelliform macular dystrophy 2-like protein 3
<b>Clone:</b>	-
<b>Applications:</b>	<b>ELISA</b> <b>WB</b>
<b>Conjugation:</b>	Non-conjugated
<b>Reactivity:</b>	Human, Mouse

### Antibody Data

---

<b>Isotype:</b>	IgG
<b>Uniprot:</b>	Q8N1M1
<b>Host Species:</b>	Rabbit
<b>Purification:</b>	Antigen Affinity purified
<b>Immunogen:</b>	Human BEST3
<b>Immunogen Species:</b>	Homo sapiens (Human)
<b>Buffer:</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Form:</b>	-

## Preparation & Storage

---

**Storage:** Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  
Store Bradford Reagent at Room Temperature for 1 Year.

<b>Recommended Dilutions:</b>	<b>Application</b>	<b>Recommended Dilution</b>
	-	-

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data

---

<b>Image</b>	<b>Description</b>
-	-