

FAM195A Antibody

PACO08081

Description

This FAM195A Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO08081
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti human PolyClonal Antibody
Synonyms:	C16orf14 antibody, c349E10.1 antibody, F195A_HUMAN antibody, FAM195A antibody, Family with sequence similarity 195 member A antibody, Protein FAM195A antibody
Clone:	-
Applications:	ELISA WB IF
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse

Antibody Data

Isotype:	IgG
Uniprot:	Q9BUT9
Host Species:	Rabbit
Purification:	Antigen Affinity purified
Immunogen:	Human C16orf14
Immunogen Species:	Homo sapiens (Human)
Buffer:	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Form:	-

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-