

LRRTM1 Antibody

PACO10295

Description

This LRRTM1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO10295
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti human PolyClonal Antibody
Synonyms:	4632401D06Rik antibody, AW125451 antibody, FLJ32082 antibody, Leucine rich repeat transmembrane neuronal 1 antibody, Leucine rich repeat transmembrane neuronal 1 protein antibody, Leucine-rich repeat transmembrane neuronal protein 1 antibody, LRRT1_HUMAN antibody, LRRTM1 antibody, Lrrtm1 leucine rich repeat transmembrane neuronal 1 antibody, OTTHUMP00000194829 antibody
Clone:	-
Applications:	ELISA WB IHC
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse, Rat

Antibody Data

Isotype:	IgG
Uniprot:	Q86UE6
Host Species:	Rabbit
Purification:	Antigen Affinity purified
Immunogen:	Human LRRTM1
Immunogen Species:	Homo sapiens (Human)
Buffer:	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Form:	-

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-