

## PRELID1 Antibody

PACO11456

### Description

---

This PRELID1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	PACO11456
<b>Contents:</b>	100µl Bradford Reagent: 1 vial (2ml)
<b>Category:</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Synonyms:</b>	25 kDa protein of relevant evolutionary and lymphoid interest antibody, CGI 106 antibody, MGC87972 antibody, mitochondrial antibody, PRELI antibody, PRELI domain containing 1 antibody, PRELI domain containing protein 1 antibody, PRELI domain containing protein 1 mitochondrial antibody, PRELI domain-containing protein 1 antibody, PRELID 1 antibody, PRELID1 antibody, PRLD1_HUMAN antibody, Protein of relevant evolutionary and lymphoid interest antibody, PX 19 antibody, PX19 antibody, Px19 like protein antibody, Px19, chicken, homolog of antibody, Px19-like protein antibody, SBBI12 antibody
<b>Clone:</b>	-
<b>Applications:</b>	<b>ELISA</b> <b>IHC</b>
<b>Conjugation:</b>	Non-conjugated
<b>Reactivity:</b>	Human, Mouse, Rat

### Antibody Data

---

<b>Isotype:</b>	IgG
<b>Uniprot:</b>	Q9Y255
<b>Host Species:</b>	Rabbit
<b>Purification:</b>	Antigen Affinity Purified
<b>Immunogen:</b>	Human PRELID1
<b>Immunogen Species:</b>	Homo sapiens (Human)

**Buffer:** PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

**Form:** -

## Preparation & Storage

---

**Storage:** Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data

---

Image	Description
-	-