

Exendin-4 Antibody

PACO19630

Description

This Exendin-4 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO19630
Contents:	50µl Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Exendin-4 antibody, Exenatide antibody
Clone:	-
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Heloderma suspectum

Antibody Data

Isotype:	IgG
Uniprot:	P26349
Host Species:	Rabbit
Purification:	Antigen affinity purification
Immunogen:	Synthetic peptide of Heloderma suspectum Exendin 4
Immunogen Species:	Heloderma suspectum (Gila monster)
Buffer:	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Form:	Liquid

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	ELISA	1:1000-1:2000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-