

CFAP410 Antibody

PACO25352

Description

This CFAP410 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO25352
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	A2 antibody, C21orf HUMF09G8.5 antibody, C21orf-HUMF09G8.5 antibody, C21orf2 antibody, Chromosome 21 open reading frame 2 antibody, CU002_HUMAN antibody, Hypothetical protein LOC755 antibody, Nuclear encoded mitochondrial protein cDNA A2 YF5 antibody, Protein C21orf2 antibody, Uncharacterized protein C21orf2 antibody, YF 5 antibody, YF5 antibody, YF5/A2 antibody
Clone:	Polyclonal
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	O43822
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Cilia- and flagella-associated protein 410 protein (1-200AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Immunohistochemistry of paraffin-embedded human gastric cancer using PACO25352 at dilution of 1:100