

IL4 Antibody

PACO27621

Description

This IL4 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO27621
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	IL4Interleukin-4 antibody, IL-4 antibody, B-cell stimulatory factor 1 antibody, BSF-1 antibody, Lymphocyte stimulatory factor 1 antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Dog

Antibody Data

Isotype:	IgG
Uniprot:	O77762
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Dog Interleukin-4 protein (25-132AA)
Immunogen Species:	Dog
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-