

POGZ Antibody

PACO29046

Description

This POGZ Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO29046
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	KIAA0461 antibody, MGC71543 antibody, Pogo transposable element with ZNF domain antibody, POGZ antibody, POGZ_HUMAN antibody, putative protein product of Nbla00003 antibody, SUHW5 antibody, suppressor of hairy wing homolog antibody, Suppressor of hairy wing homolog 5 antibody, Zinc finger protein 280E antibody, Zinc finger protein 635 antibody, ZNF280E antibody, ZNF635 antibody, ZNF635M antibody
Clone:	Polyclonal
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q7Z3K3
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Pogo transposable element with ZNF domain protein (1-300AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Immunohistochemistry of paraffin-embedded human prostate cancer using PACO29046 at dilution of 1:100

Immunohistochemistry of paraffin-embedded human kidney tissue using PACO29046 at dilution of 1:100