

lasA Antibody

PACO34142

Description

This lasA Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO34142
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	lasA antibody, PA1871Protease LasA antibody, EC 3.4.24.- antibody, Staphylytic protease antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Pseudomonas aeruginosa

Antibody Data

Isotype:	IgG
Uniprot:	P14789
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Pseudomonas aeruginosa Protease LasA protein (237-418AA)
Immunogen Species:	Pseudomonas aeruginosa (strain ATCC 15692 / DSM 22644 / CIP 104116 / JCM 14847 / LMG 12228 / 1C / PRS 101 / PAO1)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-