

SES1 Antibody

PACO34706

Description

This SES1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO34706
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	SES1 antibody, SERS antibody, YDR023W antibody, YD9813.01 antibody, Serine--tRNA ligase antibody, cytoplasmic antibody, EC 6.1.1.11 antibody, Seryl-tRNA synthetase antibody, SerRS antibody, Seryl-tRNA(Ser/Sec) synthetase antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Saccharomyces cerevisiae

Antibody Data

Isotype:	IgG
Uniprot:	P07284
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Saccharomyces cerevisiae Serine--tRNA ligase, cytoplasmic protein (1-462AA)
Immunogen Species:	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-