

## G Antibody

PACO34998

### Description

---

This G Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                      |                                        |
|----------------------|----------------------------------------|
| <b>SKU:</b>          | PACO34998                              |
| <b>Contents:</b>     | 50µg<br>Bradford Reagent: 1 vial (2ml) |
| <b>Category:</b>     | -                                      |
| <b>Synonyms:</b>     | GGlycoprotein antibody                 |
| <b>Clone:</b>        | Polyclonal                             |
| <b>Applications:</b> | <a href="#">ELISA</a>                  |
| <b>Conjugation:</b>  | Non-conjugated                         |
| <b>Reactivity:</b>   | Rabies virus                           |

### Antibody Data

---

|                           |                                                                               |
|---------------------------|-------------------------------------------------------------------------------|
| <b>Isotype:</b>           | IgG                                                                           |
| <b>Uniprot:</b>           | A3RM22                                                                        |
| <b>Host Species:</b>      | Rabbit                                                                        |
| <b>Purification:</b>      | >95%, Protein G purified                                                      |
| <b>Immunogen:</b>         | Recombinant Rabies virus Glycoprotein protein (20-459AA)                      |
| <b>Immunogen Species:</b> | Rabies virus                                                                  |
| <b>Buffer:</b>            | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| <b>Form:</b>              | Liquid                                                                        |

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

## Preparation & Storage

---

**Storage:** Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  
Store Bradford Reagent at Room Temperature for 1 Year.

| Recommended Dilutions: | Application | Recommended Dilution |
|------------------------|-------------|----------------------|
|                        | -           | -                    |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data

---

| Image | Description |
|-------|-------------|
| -     | -           |