

lepB Antibody

PACO35226

Description

This lepB Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO35226
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	lepB antibody, Rv2903c antibody, MTCY274.34c antibody, Signal peptidase I antibody, SPase I antibody, EC 3.4.21.89 antibody, Leader peptidase I antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Mycobacterium tuberculosis

Antibody Data

Isotype:	IgG
Uniprot:	P9WKA1
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Mycobacterium tuberculosis Signal peptidase I protein (88-294AA)
Immunogen Species:	Mycobacterium tuberculosis
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-