

TREML2 Antibody

PACO36538

Description

This TREML2 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO36538
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	TREML2 antibody, C6orf76 antibody, TLT2 antibody, UNQ6268/PRO20473 antibody, Trem-like transcript 2 protein antibody, TLT-2 antibody, Triggering receptor expressed on myeloid cells-like protein 2 antibody
Clone:	Polyclonal
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q5T2D2
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Trem-like transcript 2 protein (19-268AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Immunohistochemistry of paraffin-embedded human colon cancer using PACO36538 at dilution of 1:100

Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO36538 at dilution of 1:100