

ARHGAP27 Antibody

PACO36718

Description

This ARHGAP27 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO36718
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	ARHGAP27 antibody, CAMGAP1 antibody, CIN85-associated multi-domain-containing Rho GTPase-activating protein 1 antibody, PP905 antibody, RHG27_HUMAN antibody, Rho GTPase-activating protein 27 antibody, Rho-type GTPase-activating protein 27 antibody, SH3 domain containing 20 antibody, SH3 domain containing protein 20 antibody, SH3D20 antibody, SH3P20 antibody
Clone:	Polyclonal
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q6ZUM4
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Rho GTPase-activating protein 27 protein (1-300AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol


Validation Data

Image



Description

Immunohistochemistry of paraffin-embedded human lung cancer using PACO36718 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic cancer using PACO36718 at dilution of 1:100