

MACROD1 Antibody

PACO40286

Description

This MACROD1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO40286
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	MACROD1 antibody, LRP16 antibody, ADP-ribose glycohydrolase MACROD1 antibody, MACRO domain-containing protein 1 antibody, O-acetyl-ADP-ribose deacetylase MACROD1 antibody, EC 3.1.1.106 antibody, Protein LRP16 antibody, [Protein ADP-ribosylaspartate] hydrolase MACROD1 antibody, EC 3.2.2.- antibody, [Protein ADP-ribosylglutamate] hydrolase MACROD1 antibody, EC 3.2.2.- antibody
Clone:	Polyclonal
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q9BQ69
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human O-acetyl-ADP-ribose deacetylase MACROD1 protein (83-325AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Immunohistochemistry of paraffin-embedded human colon cancer using PACO40286 at dilution of 1:100