

CHMP4C Antibody

PACO46994

Description

This CHMP4C Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO46994
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Charged multivesicular body protein 4c antibody, CHM4C_HUMAN antibody, CHMP4c antibody, Chromatin-modifying protein 4c antibody, hSnf7-3 antibody, hVps32-3 antibody, SHAX3 antibody, SNF7 homolog associated with Alix 3 antibody, SNF7-3 antibody, Vacuolar protein sorting-associated protein 32-3 antibody, Vps32-3 antibody
Clone:	Polyclonal
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q96CF2
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Charged multivesicular body protein 4c protein (1-233AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Immunohistochemistry of paraffin-embedded human colon cancer using PACO46994 at dilution of 1:100