

HA Antibody

PACO50462

Description

This HA Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO50462
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	HA antibody, Hemagglutinin [Cleaved into: Hemagglutinin HA1 chain, Hemagglutinin HA2 chain] antibody
Clone:	Polyclonal
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Influenza A virus

Antibody Data

Isotype:	IgG
Uniprot:	Q0HD60
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Influenza A virus Hemagglutinin protein (18-343AA)
Immunogen Species:	Influenza A virus
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:500-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western Blot Positive WB detected in Recombinant protein All lanes: HA antibody at 3µg/ml Secondary Goat polyclonal to rabbit IgG at 1/50000 dilution predicted band size: 64 kDa observed band size: 64 kDa