

Anti-Ipilimumab ADA ELISA Kit

SKU: GXEL00148 | Assay: Quantitative ELISA (Colorimetric)

The Anti-Ipilimumab ADA ELISA Kit (SKU: GXEL00148) is a high-sensitivity quantitative ELISA for the quantitative measurement of Anti-Ipilimumab in Plasma, Serum. The assay covers a detection range of 156.25 - 10, 000 ng/mL and is read colorimetrically at 450 nm.

Specifications

Product Name	Anti-Ipilimumab ADA ELISA Kit
SKU	GXEL00148
Accession	P16410
Alternative Names	BMS-734016, MDX-010, MDX-CTLA-4, 477202-00-9
Assay Type	Quantitative
Detection Method	Colorimetric
Detection Range	156.25 - 10, 000 ng/mL
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00148.pdf> | Version 1.0, 09 July 2026