

Bovine LIF Recombinant Protein

Description

This high-purity Bovine LIF Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

| | |
|-----------------------------|--|
| SKU: | KPRT0274 |
| Contents: | Protein: 1 vial (5µg) Bradford Reagent: 1 vial (2ml) |
| Source: | Yeast |
| Application: | Cell Culture, ELISA, ELISpot, WB |
| Molecular Weight: | 19.8 kDa |
| Purity: | >98% as visualized by SDS-PAGE analysis |
| Amino Acid Sequence: | SPLPITPVNA TCATRHPCPS NLMNQIRNQL GQLNSSANSI FILYYTAQGE PFPNNLDKLC SPNVTDFFPPF HANATEKARL VELYRIIAYL GASLGNITRD QKVLNPNYAHG LHSKLSTTAD VLRGLLSNVL CRLCSKYHVS HVDVITYGPD SGKDVFQKKK LGCQLLGKYK QVIAVLAQAF |
| Homology: | Bos taurus (cattle) |

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.