

Turkey IL-8 Recombinant Protein

Description

This high-purity Turkey IL-8 Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT0306
Contents:	Protein: 1 vial (5µg / 25µg) Bradford Reagent: 1 vial (2ml)
Source:	Yeast
Application:	Cell Culture, ELISA, ELISpot, WB
Molecular Weight:	17.3 kDa
Purity:	>98% as visualized by SDS-PAGE analysis
Amino Acid Sequence:	QIGHPSPPSE KKEKKAHL TGKLNRSIP LEWEDTYGIA LVSGVKYKKG SLVINETGLY FVYSKVYFRG QSCNNQPLSH KVYSRNFYYP QDMVLMGKM MNYCTSGKMW ARSSYLGAVF NLTSADHLYV NVSELSLVSF EESKTFFGLY KL
Homology:	Bison bison bison (American buffalo), Bos indicus (zebu), Bos taurus (cattle), Bubalus bubalis (water buffalo)

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.