

## Bovine IL-1 Receptor Antagonist Recombinant Protein

### Description

---

This high-purity Bovine IL-1 Receptor Antagonist Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Protein Information

---

<b>SKU:</b>	KPRT0351		
<b>Contents:</b>	Protein: 1 vial (5µg / 25µg / 100µg / 500µg(5X100µg)) Bradford Reagent: 1 vial (2ml)		
<b>Source:</b>	Yeast		
<b>Application:</b>	Cell Culture, ELISA, ELISpot, WB		
<b>Molecular Weight:</b>	8.6 kDa		
<b>Purity:</b>	>98% as visualized by SDS-PAGE analysis		
<b>Amino Acid Sequence:</b>	QQHGVNKCNL	FCNKMTQKIP	ESRLVGYQRN RESCNDGAVI
	LKTVKGKSFC	ADPKEEWVQK	AMKHLDHKAS VSQKSG
<b>Homology:</b>	Bos indicus (zebu), Bos taurus (cattle)		

## Preparation & Storage

---

**Stability & Storage:** Store Protein at -20°C for 1 Year.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.