

## Porcine SCF Recombinant Protein

### Description

---

This high-purity Porcine SCF Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Protein Information

---

<b>SKU:</b>	KPRT0391
<b>Contents:</b>	Protein: 1 vial (5µg) Bradford Reagent: 1 vial (2ml)
<b>Source:</b>	Yeast
<b>Application:</b>	Cell Culture, ELISA, ELISpot, WB
<b>Molecular Weight:</b>	18.7 kDa
<b>Purity:</b>	>98% as visualized by SDS-PAGE analysis
<b>Amino Acid Sequence:</b>	QGICRNRVTD DVKDVTKLVA NLPKDYKITL KYVPGMDVLP SHCWISEMVE QLSVSLTDLL DKFSNISEGL SNYSIIDKLV KIVDDLVECM EEHSFENVRK SSKSPEPRLF TPEKFFGIFN RSIDAFKDLE MVAPKTSECV ISSTLTPEKD SRVSVTKPFM LPPVAA
<b>Homology:</b>	Sus scrofa (pig)

## Preparation & Storage

---

**Stability & Storage:** Store Protein at -20°C for 1 Year.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.