

Cynomolgus Monkey CXCL4 Recombinant Protein

Description

This high-purity Cynomolgus Monkey CXCL4 Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT0589
Contents:	Protein: 1 vial (5µg) Bradford Reagent: 1 vial (2ml)
Source:	Yeast
Application:	Cell Culture, ELISA, ELISpot, WB
Molecular Weight:	7.8 kDa
Purity:	>98% as visualized by SDS-PAGE analysis
Amino Acid Sequence:	AEAEEDGDLQ CLCVKTTSQV RPRHITSLEV IKAGPHCPTA QLIATLKNGR KICLDLQAPL YKKIVKKLLE S
Homology:	Cercocebus atys (sooty mangabey), Chlorocebus sabaeus (green monkey), Macaca fascicularis (crab-eating macaque), Macaca mulatta (Rhesus monkey), Macaca nemestrina (pig-tailed macaque), Mandrillus leucophaeus (drill), Piliocolobus tephrosceles (Ugandan red Colobus), Rhinopithecus bieti (black snub-nosed monkey), Rhinopithecus roxellana (golden snub-nosed monkey), Theropithecus gelada (gelada)

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.