

## Cynomolgus Monkey GM-CSF Recombinant Protein

### Description

---

This high-purity Cynomolgus Monkey GM-CSF Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Protein Information

---

<b>SKU:</b>	KPRT1185
<b>Contents:</b>	Protein: 1 vial (25µg, 100µg, 500µg(5X100µg)) Bradford Reagent: 1 vial (2ml)
<b>Source:</b>	Yeast
<b>Application:</b>	Cell Culture, ELISA, ELISpot, WB
<b>Molecular Weight:</b>	14.4 kda
<b>Purity:</b>	>98% as visualized by SDS-PAGE analysis
<b>Amino Acid Sequence:</b>	APARSPSPGT QPWEHVNAIQ EARRLLNLSR DTAAEMNKTV EVSSEMF DLQ EPSCLQTRLE LYKQGLQGS L TKLKGPLTMM ASHYKQHCPP TPETSCATQI ITFQSFKENL KDFLLVIPFD CWEPVQE
<b>Homology:</b>	Macaca fascicularis (crab-eating macaque), Macaca mulatta (Rhesus monkey)

## Preparation & Storage

---

**Stability & Storage:** Store Protein at -20°C for 1 Year.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.