

Rabbit TGF-alpha Recombinant Protein

Description

This high-purity Rabbit TGF-alpha Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT1648		
Contents:	Protein: 1 vial (25µg, 100µg, 500µg(5X100µg)) Bradford Reagent: 1 vial (2ml)		
Source:	Yeast		
Application:	Cell Culture, ELISA, ELISpot, WB		
Molecular Weight:	5.6 kDa		
Purity:	>98% as visualized by SDS-PAGE analysis		
Amino Acid Sequence:	VVSHFNQCPD	SHTQFCFHGT	CRFLVQEDKP ACVCHSGYVG ARCEHADLLA
Homology:	Octodon degus (degu), Oryctolagus cuniculus (rabbit)		

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.