

Chicken/Mallard FGF Basic Recombinant Protein

Description

This high-purity Chicken/Mallard FGF Basic Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT2122
Contents:	Protein: 1 vial (5µg, 25µg, 100µg, 500µg(5X100µg)) Bradford Reagent: 1 vial (2ml)
Source:	Yeast
Application:	Cell Culture, ELISA, ELISpot, WB
Molecular Weight:	16.3 kDa
Purity:	>98% as visualized by SDS-PAGE analysis
Amino Acid Sequence:	STPEDGGSGG FTPGGFKDPK RLYCKNGGFF LRIRSDGGVD GIREKSDAHI KLQIQATSVG EVVIKGVCAN RYLAMNRDGR LFGVRRATDE CYFLERLESN NYNTYRSRKY PGMVYALKRT GQYKSGSKTG PGQKAILFLP MSAKC
Homology:	Anas platyrhynchos (mallard), Anser cygnoides (Swan goose), Aythya fuligula (tufted duck), Centrocercus urophasianus (Greater sage-grouse), Coturnix japonica (Japanese quail), Cygnus atratus (black swan), Cygnus olor (mute swan), Gallus gallus (chicken), Lagopus leucura (white-tailed ptarmigan), Lagopus muta (rock ptarmigan), Oxyura jamaicensis (ruddy duck), Phasianus colchicus (Ring-necked pheasant)

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.