

RPCB1102

---

## Product Information

<b>Product SKU:</b>	RPCB1102	<b>Gene ID:</b>	928966	<b>Size:</b>	20µg
<b>Tag:</b>	C-His	<b>Reactivity:</b>	Monkeypox virus		

---

## Additional Information

<b>Expression Host:</b>	E. coli	<b>Swissprot:</b>	Q77HM6
<b>Purity:</b>	> 85% by SDS-PAGE.		

---

## Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Monkeypox virus A29 Protein , tested reactivity in E. coli and has been validated in SDS-PAGE. 100% guaranteed.

**Endotoxin:** <0.1 EU/µg

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

**Storage:** Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.