

RPCB1197

Product Information

Product SKU:	RPCB1197	Gene ID:	6356	Size:	10µg
Tag:	C-hFc	Reactivity:	Human		

Additional Information

Expression Host:	HEK293 cells	Swissprot:	P51671
Purity:	> 95% by SDS-PAGE.		

Protein Information

Background: In response to the presence of allergens, this protein directly promotes the accumulation of eosinophils, a prominent feature of allergic inflammatory reactions

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human CCL11/Eotaxin Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: <0.1EU/µg

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.