

RPCB1536

Product Information

Product SKU: RPCB1536 **Gene ID:** 54790 **Size:** 10µg
Tag: Please contact us for more information. **Reactivity:** Human

Additional Information

Expression Host: HEK293 cells **Swissprot:** Q6N021-1
Purity: > 90% by SDS-PAGE.

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human TET2 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: Please contact us for more information.

Formulation: Supplied as a 0.22 µm filtered solution of 50mM HEPES,300mM NaCl,1mM TCEP,10% Glycerol,pH 8.0.

Storage: Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.