

RPCB1883

Product Information

Product SKU:	RPCB1883	Gene ID:	16854	Size:	100µg
Tag:	C-His	Reactivity:	Mouse		

Additional Information

Expression Host:	HEK293 cells	Swissprot:	Q8C253
Purity:	> 95% by SDS-PAGE.		

Protein Information

Background: Galectin-3 (Gal-3) is a 30 kDa β -galactose, highly conserved and widely distributed intracellularly and extracellularly. Gal-3 has been demonstrated in recent years to be a novel inflammatory factor participating in the process of intravascular inflammation, lipid endocytosis, macrophage activation, cellular proliferation, monocyte chemotaxis, and cell adhesion.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Mouse Galectin-3/LGALS3 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 1 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.